


NSE directly measures the intermediate scattering function $I(Q, t)$

The NIST-NSE spectrometer is best used for measuring coherent diffusive or dispersionless excitations at long time in the range of 0.01 to 100 ns (several neV to several hundred μ eV) and 0.02 to 1.5 \AA^{-1} of the length scale.

Problems in the dynamical range are

1. Polymers
Observation and quantitative description of the crossover in dynamics from local segmental diffusion to time-dependent behavior governed by entanglements occurring over longer length scale
2. Glassy dynamics
Identification in polymer glasses of the intra- and inter-molecular dynamics responsible for the α and β_{slow} relaxation
3. Biological model systems
Quantitative description of the effect of interlayer coupling in the extended diffusive mode of lipid bilayers
4. Proteins
Intra-molecular diffusion in e.g., pig immunoglobulin G
5. Slow magnetic spin dynamics
Spin dynamics in spin-glasses, frustrated magnets and magnetoresistive materials for example

Other forms of dynamical behavior (in the same dynamical range) are also accessible, with some more difficulty, such as incoherent dynamics.




NSE covers the dynamical range in between conventional inelastic neutron scattering and dynamic light scattering techniques

NSE principle


NSE technique has the highest energy resolution among the inelastic neutron scattering techniques. In the case of the conventional inelastic neutron scattering techniques, the velocities of the incident and the scattered neutrons have to be determined in order to know the energy transfer of neutrons at the sample position. Comparing to the conventional techniques, NSE technique uses the Larmor precession of the neutron spin in the magnetic field as a measure of the energy transfer of neutrons. This unique feature allows to decouple the energy resolution of neutron itself from the dynamical modes in the measuring sample. This is the reason why the highest energy resolution has achieved by this technique without losing the neutron intensity.

Comparison between NSE and conventional inelastic neutron scattering techniques

Neutron spin echo technique



Conventional inelastic neutron scattering techniques



Characteristics of NG5-NSE spectrometer

wavelength, λ	5 -- 12 \AA
scattering angle, 2θ	0 -- 105 °
maximum field integral, $\int B \cdot dl$	0.438 Tm
field inhomogeneity, $\Delta J/J$	$\sim 10^{-6}$
Fourier time, t	0.005 -- 0.1 ns (by shorty operation) 0.05 -- 100 ns (by normal operation)
Momentum transfer, q	$q_{\min} \sim 0.02 \text{\AA}^{-1}$ (by 11 \AA) $q_{\max} \sim 1.66 \text{\AA}^{-1}$ (by 6 \AA)
beam size at sample	30 x 30 mm ²
polarizer	Mezei V cavity (transmission polarizer)
analyzer	3q _c Fe/Si supermirror Pi192 Co/Ti supermirrors
flipping ratio, Up/down	9 for 6 \AA 8 for 8 \AA 7 for 11 \AA

