



# Evaluation of thiol acetates ( $\text{CH}_3\text{OSR}$ ) and disulfides ( $\text{RSSR}$ ) as lipidic anchors in tethered bilayer membranes (tBLMs)

Perri Weinstein

Advisors: Dr. David Vanderah and Dr. Frank  
Heinrich

# Integral Membrane Proteins (IMPs)


- Cell membranes composed of phospholipid bilayer
  - Hydrophilic and hydrophobic
- Proteins carry out a variety of functions
  - Channels
  - Receptors
  - Enzymes
- Structure and function studies



# Why we are doing this work

- Native bilayers are complex
  - Many different types of IMPs and lipids
- Structure and function studies are difficult with IMPs
- Tethered Bilayer Lipid Membranes (tBLMs) create a model for studying IMPs
  - Environment stable for days
  - Surface measurements


# tBLMs



# tBLM Preparation – Rapid Solvent Exchange Method

## Step 1: Mixed SAM

Incubation with tether & small molecule




Ng 2008 SURF talk

- Gold substrates sit in solution for 1-5 days

## Step 2: Create tBLM

Rapid Solvent Exchange



- Floating bilayer to replicate cell membrane fluidity



# The Key Players

- **WC<sub>14</sub>**

- **HS**(CH<sub>2</sub>CH<sub>2</sub>O)<sub>6</sub>CH<sub>2</sub>CH(OR)CH<sub>2</sub>(OR)
  - where R=C<sub>14</sub>H<sub>29</sub>

- [β-mercaptoethanol (**βME**)]

- **HS**CH<sub>2</sub>CH<sub>2</sub>OH


- Thiols lack long term stability!!!!!! Looking for alternatives

- **Thiol acetates**

- **H<sub>3</sub>COS**(CH<sub>2</sub>CH<sub>2</sub>O)<sub>6</sub>CH<sub>2</sub>CH(OR)CH<sub>2</sub>(OR)
- **H<sub>3</sub>COS**CH<sub>2</sub>CH<sub>2</sub>OCH<sub>3</sub> (βME methyl ether = βMEOMe)


- **Disulfides**

- R'**SS**R', where R' = (CH<sub>2</sub>CH<sub>2</sub>O)<sub>6</sub>CH<sub>2</sub>CH(OR)CH<sub>2</sub>(OR)
- HOCH<sub>2</sub>CH<sub>2</sub>**SS**CH<sub>2</sub>CH<sub>2</sub>OH



# Surface Measurements


- Ellipsometry
  - An optical metrology that allows the determination of thin film thicknesses
- Contact Angle
  - A technique that allows the determination of the surface energy (hydrophobicity/hydrophilicity)
- Neutron Reflectivity
  - A scattering technique that allows the determination of the structure of films
- Electrochemical Impedance Spectroscopy
  - An electrochemical technique that measures the resistance and capacitance of films as a function of frequency to describe the bilayer and its imperfections



# Ellipsometry Data: thiol acetates

- Refresher: determination of thin film thicknesses
- Graph of 1 of 3 data sets
- SAM growth as a function of increasing tether concentration and time

**WC<sub>14</sub> Thiol acetate &  $\beta$ MEOMe Solutions at 0.2 mM**






# Ellipsometry Data: $R'SSR'/[\beta MES]_2$

- Graph of 1 of 2 data sets
- SAM growth as a function of increasing tether concentration and time
- Anomalies could be due to large molecules laying down on the surface → “blanket theory”


**WC<sub>14</sub> Disulfide &  $\beta$ ME Disulfide Solutions at 0.1 mM**



# Contact Angle Data: thiol acetates

- Refresher: determination of the hydrophobicity/hydrophilicity of the surface
- Contact Angle less than WC<sub>14</sub> thiol → surface is less ordered


**Contact Angle Measurements of Thiol Acetates**



# Contact Angle Data: disulfides


- Same results for disulfides as thiol acetates: not as ordered as thiols

**Contact Angle Measurements of Disulfides**




# Neutron Reflectometry

- Refresher: determination of structure of films based on scattering of a reflected neutron beam
- Using a known scattering length of each nuclei and a known number of atoms
  - calculate material density (SLD)
- Sample: 70:30 tether composition  $\rightarrow$  50% less SAM seen than WC<sub>14</sub>,
  - NOT well ordered
- Only incubated for 24 hours
  - Would have told us about the sub-membrane space if timing was not an issue




# Electrochemical Impedance Spectroscopy (EIS)

## WC<sub>14</sub> Capacitance



## WC<sub>14</sub>SAc Capacitance (100%)



- Refresher: measures the resistance and capacitance of films as a function of frequency
- Performed on 4 samples of WC<sub>14</sub>SAc:βMEOME: (A) 20: 80 (B) 30:70 (C) 40:60, (D) 100:0 (3 days)
  - Lower concentrations: bilayers did not form consistently
  - 100:0 SAM is reasonable and indicates good coverage, 100:0 BLM (bottom curve)
  - Reasonable, but higher capacitance than obtain with thiols → leaky (less ordered, more defects)
- Overall, there's variability however the data indicates the formation of SAMs and tBLMs with thiol acetates

# Conclusions and Upcoming Research


- Can thiol acetates and/or disulfides replace thiols in tBLMs?
- Maybe...
  - Less ordered
    - Lower contact angles
    - Increased capacitance
- So we are close but no cigar...yet
  - Collect more data on thiol acetates and disulfides
  - Especially with disulfides
    - Thicknesses very comparable to WC14 at high tether concentrations
- However, demonstrated the formation of a tBLM formation with thiol acetates
  - Look at long term stability



# Acknowledgments

- David Vanderah
- Frank Heinrich
- Sandra Gutierrez
- Marlon Walker
- David Hoogerheide
- Duncan McGillivray
- John Marino

thank  
you!



# Questions?

