

# MACS Concept and Project Status

Collin Broholm

Johns Hopkins University and NIST Center for Neutron Research



- Making best use of CW source
- MACS-imizing incident flux
- MACS-imizing detection efficiency
- From concept to reality

# Goals in Neutron Spectroscopy


---

- **A technique with a future**
  - Unique information about dynamic correlations
  - Model independent access to interaction strength
  - Access microscopic structure of dynamic systems
- **Limited scope on current instruments**
  - Need cm<sup>3</sup> sized crystals
  - Need weeks of beam time
- **Increased sensitivity will broaden impact**
  - Smaller samples earlier in new materials cycle
  - Parametric studies as in diffraction
  - Comprehensive surveys for tests of theory


# The Spallation Neutron Source





# MAPS Spectrometer at ISIS



# Spectroscopy from MAPS



$KCuF_3$   
*A. Tennant et al.*



$Y_{2-x}Ca_xBaNiO_5$   
*Y. Chen et al*

# Advantages of TAS like instrumentation

---

- Use Bragg optics to focus incident beam
- Use Bragg reflection to polarize beam
- **Select probed energy range  
at fixed resolution and sensitivity**
- Flexibility in resolution for given "kinematic footprint"

**TAS sensitivity scales with CW flux so are best implemented at Reactor based neutron source`**

# Why cold neutrons and double focusing

- Q and E resolved spectroscopy requires

$$\Delta E \approx 0.1 J \quad \Delta Q \approx 0.1 a^{-1}$$

- **Energy scale  $J$  varies more than length scale  $a$**

Lattice	Compound	J (meV)	a (Å)
3D Cubic S=5/2	La <sub>0.7</sub> Pb <sub>0.3</sub> MnO <sub>3</sub>	8.8	3.9
2D Square S=1/2	La <sub>2</sub> CuO <sub>4</sub>	132	5.4
2D Kagomé S=3/2	KCr <sub>3</sub> (OH) <sub>6</sub> (SO <sub>4</sub> ) <sub>2</sub>	1.2	3.7
2D Kagomé S=3/2	SrCr <sub>9</sub> Ga <sub>3</sub> O <sub>19</sub>	10	2.9
1D S=1/2 chain	Cu(C <sub>6</sub> D <sub>5</sub> COO) <sub>2</sub> ·3D <sub>2</sub> O	1.5	3.2
1D S=1/2 chain	KCuF <sub>3</sub>	35	3.9
1D S=1 chain	NENP	4.1	5.2
1D S=1 chain	AgVP <sub>2</sub> S <sub>6</sub>	58	2.9

- To probe a range of materials must vary


$$10^{-2} \text{ meV} < \Delta E < 10 \text{ meV} \quad \text{keeping } \Delta Q \approx 0.05 \text{ \AA}^{-1}$$

- To probe hard matter with low energy scales

- Reduce  $E_i$ . Cold source provides the flux

- Increase angular divergence before and after sample


# NCNR Liquid Hydrogen cold source




New cold source to be installed in 2001 will double flux




# Bragg focusing






# Projected performance analytical approximation



# Monte Carlo Simulation of MACS




# Increasing Efficiency of the detection system

---

- Maximize solid angle of detection
- Maintain the virtues of TAS
  - Ability to select single energy
  - Ability to vary resolution at fixed  $E_f$
  - High signal to noise ratio
- Detect neutrons "not at  $E_f$ "
- Need a practical solution

# Double bounce analyzer array




Energy integrating  
detector


Spectroscopic detector

# Constant energy transfer slice

$$E_f = 3.7 \text{ meV} \quad \hbar\omega = 1 \text{ meV}$$




# Assembling slices to probe Q-E volume





# Fixed vertical focusing of analyzers



Double analyzer is "compound lens" → efficient vertical focusing

# The history of MACS

---

- 1993** Discussions about the possibility of a "sub-thermal" TAS on NGO
- 1994** Analytical calculations show efficacy of double focusing at NGO
- 1995** Initiate JHU/NIST project to develop conceptual design
- 1998** Top level specification for monochromator completed
- 1998** JHU/NIST project starts to develop Doubly Focusing Monochromator
- 2000** Christoph Brocker starts engineering design

# Scientific Program for NGO spectrometer

---

- **Expanding the scope for Inelastic scattering from crystals:**
  - 0.5 mm<sup>3</sup> samples
  - impurities at the 1% level
  - complete surveys to reveal spin-wave-conduction electron interactions
  - extreme environments: pressure and fields to tune correlated systems
- **Probing short range order**
  - solid ionic conductors, spin glasses, quasi-crystals, ferroelectrics, charge and spin polarons, quantum magnets, frustrated magnets.
- **Weak broken symmetry phases**
  - Incommensurate charge, lattice, and spin order in correlated electron systems
- **Excitations in artificially structured solids**
  - Spin waves in magnetic super-lattices
  - magnetic fluctuations in nano-structured materials

# MACS development team

## Overall instrument design

Paul Brand	NIST
Christoph Brocker	NIST
Collin Broholm	JHU
Jeff Lynn	NIST
Mike Rowe	NIST
Jack Rush	NIST
Yiming Qiu	JHU

## Shielding Calculations

Charles M. Eisenhauer NIST

## Focusing Monochromator

Steve Conard	JHU
Joe Orndorff	JHU
Tim Reeves	JHU
Gregg Scharfstein	JHU
Stephen Smee	UMD



# Conclusions

---

- Large solid angle access to NCNR cold source enables a unique cold neutron spectrometer.
- MACS employs Bragg optics to focus the beam and provide  $> 10^8$  n/cm<sup>2</sup>/s on the sample at 0.2 meV resolution.
- The MACS detector concept offers the versatility, resolution, and low background of a TAS with 20 times greater data rate.
- The highly efficient constant E mode of MACS will provide a unique capability for probing the structure of fluctuating systems.
- MACS will be complementary to the DCS and future SNS TOF spectrometers because of its unique data collection protocol
- NSF partnership is needed to realize the potential and make the unique capabilities available to widest possible cross section of the US scientific community